



Universal Service Dental Lab

Rx REMOVABLE APPLIANCE

132 W Verdugo Ave,
Burbank, CA 91502

• 888-805-0808 • 747-268-0808

smilesspecialist@usdentallab.net

Dr. Name _____ Phone # _____

Acct. # _____ Patient Name _____
First Last

Address/Email _____ Deliver by 5 p.m. on _____

Enclosed with Case: ☐ Impressions ☐ Models ☐ Bite ☐ Photos ☐ Other: _____

FINAL SHADE



Indicate Shade Here

PRESENT TOOTH OR STUMP SHADE



Indicate Shade Here

OCCUSAL STAINING

☐ None ☐ Light* ☐ Medium ☐ Dark

PONTIC DESIGN



☐ Mod RidgeLap ☐ RidgeLap ☐ Ovate ☐ Bullet ☐ Hygienic

MARGIN AND METAL DESIGN



☐ Labial Butt ☐ 360° Butt ☐ Junction ☐ Junction




NIGHTGUARDS/RETAINERS

☐ Upper Arch ☐ Lower Arch

- ☐ Comfort3D (3D-printed, hard)
☐ Comfort H/S (hard)
☐ CLEARsprint (self-adjusting, hard)
☐ Soft nightguard
☐ Playsafe mouthguards

- ☐ quantity of 2 ☐ quantity of 4 ☐ quantity of 6
☐ Clear Retainer ☐ Hawley Retainer
☐ Essix Retainer (1-tooth) ☐ Bleaching Tray



DENTURES/FLIPPERS/FLEXIBLE PARTIALS

Denture

- ☐ Handcrafted
☐ Digital (3D-printed)
☐ Immediate
 ☐ 3D-printed
 ☐ Handcrafted
☐ Copy (3D-printed)

Partial

- ☐ Flipper
☐ tcs

Select Phase

- ☐ Custom tray
☐ Bite rim
☐ Setup try-in
☐ Finish

☐ Digital Teeth Shade _____ Mould _____

☐ Premium Brand (Yamahachi)

Shade _____ Mould _____

☐ Other Premium Brand (extra charge)

Shade _____ Brand _____ Mould _____

Tooth Setup

- ☐ Ideal ☐ Characterized ☐ Study model
☐ Masculine ☐ Feminine ☐ Age _____

Gingival Shade

- ☐ Std. G1
☐ Med. G3
☐ Dark G4

Flexible Partial Shade

- ☐ Lt. Pink ☐ Lt./Dark Pink
☐ Std. ☐ Dark Pink

*Standard unless specified otherwise.

SIMPLY NATURAL METAL PARTIALS

Metal frame with acrylic and Yamahachi teeth

Phase

- ☐ Metal frame try-in
☐ Printed frame try-in
☐ Frame w/occlus. rim
☐ Frame w/setup try-in
☐ Finish
☐ Scan/Save File
(extra charge applies)

☐ CAD/CAM-milled acetal partial

Color: _____

COMBINATION CROWNS & PARTIALS

- ☐ Fabricate RPD to fit restoration
☐ Future RPD
 ☐ Valplast
 ☐ Attachments
☐ Oxidian Non-Precious
☐ Oxidian Noble
 ☐ UniZir Full-Strength* (> 1,000 MPa)
 ☐ UniZir Esthetic (870 MPa)

SNORING/SLEEP APNEA APPLIANCES

(Upper and lower models with protrusive bite required)

- ☐ Silent Nite Sleep Appliance* ☐ TAP 3 TL
☐ Silent Nite with Universal Hinge ☐ dreamTAP
☐ OASYS Hinge Appliance ☐ flexTAP
☐ EMA

2020 ©. Property of Universal Service Dental.

TERMS AND CONDITIONS

All Restorations Made in the USA



TERMS: Cost of collection of any account will be paid by the customer. All accounts are payable within 15 days of statement date. **Accounts not paid within the stated terms will be subject to COD status and a late charge of 3 percent of the unpaid balance.** \$59 NSF fee will apply for returned checks. Prices subject to change without notice. Rx must be enclosed with original case submission.

NO-FAULT REMAKE POLICY: USD Lab is pleased to process all remakes or adjustments at no additional charge if requested within the warranty period and accompanied by the return of the original appliance.

LIMITED WARRANTY/LIMITATION OF LIABILITY. For warranty terms and conditions and limitation of liability, visit www.USDentallab.net

PREPARATION GUIDELINES

Rest Preparations

Occlusal Rest



1/3 width of faciolingual,
1/2 width of cusp tips




Rest depth at least 1 mm

Channel Rest



Inverted V Rest




Guide Plane



Retentive Prep



Interproximal Preparation




Occlusal rests only



Rests with
buccal/lingual access


Impressions/Model

Maxillary Arch



1. Incisive papilla
2. Labial frenum
3. Buccal frenum
4. Maxillary tuberosity
5. Hamular notch

Mandibular Arch



6. Labial frenum
7. Buccal frenum
8. Retromolar pad
9. Mylohyoid ridge

Recommended Impression Materials:

- Vinyl polysiloxane (Capture®, Imprint™, Take 1™, Aquasil, Splash!®, etc.) – *light, regular or monophasic viscosities recommended*
- Polyether (Impregum™, Permadyne™)